

# Universal Protein Isolation Kit – Peptide (UPECK-Pep) Guide

For extraction, cleanup, and digestion of protein samples

## Contents:

<b><u>Kit component</u></b>	<b><u>Amount</u></b>	<b><u>Catalog number</u></b>
ProMTag	25 µL	Pep-3001
mP-Trypsin	175 µL	Pep-3002
Lysis buffer (LB)	25 mL	Pep-3003
DTT	1	Pep-3004
IAA	8	Pep-3005
Wash buffer 1 (WB1)	8 mL	Pep-3006
Wash buffer 2 (WB2)	8 mL	Pep-3007
Wash buffer 3 (WB3)	8 mL	Pep-3008
Elution buffer (EB)	1 mL	Pep-3009
ProMTag capture tubes	8	Pep-3010
2 mL waste collection tubes	8	-
1.5 mL protein low bind tubes	16	-

## Storage:

Store the entire kit at 4°C, with the exception of the Lysis Buffer which should be kept at room temperature. Make note of the expiration date on the box. Past this date, certain elements of the box may begin to lose their efficacy. If you can't find the expiration date, email us at [info@impactproteomics.com](mailto:info@impactproteomics.com) with the lot number and we will find it for you.

## Safety:

Always protect yourself appropriately when working with chemicals. This includes, but is not limited to, an appropriate lab coat, disposable gloves, and protective eye goggles. For more information, please read the included Safety Data Sheets about the included reagents. These are also available online at [www.impactproteomics.com/X](http://www.impactproteomics.com/X)

Note: The ProMTag, WB1, WB2, and WB3 contain various amounts of acetonitrile. Please dispose of appropriately and avoid open flames.

Note: EB contains formic acid. Avoid contact with skin and eyes. If skin contact occurs, remove contaminated garments and rinse skin thoroughly with water. If eye contact occurs, remove contact lenses if applicable and rinse with plenty of water.

## Equipment and reagents you will need before you start:

- Protein sample, cell lysate, or protein source
- Pipettes and pipette tips
- Benchtop centrifuge (mini or full size)
- Sample rotator (rotisserie or carousel)
- Heating block
- Boiling water bath (if using our lysis method)
- Vortex

## Cell lysis and preparation of the biological sample for UPECK processing

No matter what biological source our kits will work to separate proteins from any other undesirable molecules. For the best results, cell lysis must be as thorough as possible. We have included our favorite lysis buffer for you to use that is compatible with most tissue lysis methods, but it is certainly not the only one that is compatible.

Feel free to use your own lysis buffer **as long as it does not contain TRIS (or any other buffer with primary amines) and is ~ pH 8.0**. If your lysis technique uses TRIS, we recommend switching to 100 mM HEPES, pH 8.0 (the concentration matters so make sure it's 100 mM).

If you are not working with a strong denaturant such as high concentration SDS, we **highly** recommend using a protease inhibitor in your lysis buffer to prevent protein degradation.

## Other notes to consider before you begin

- 1) All centrifugation steps may be performed on a benchtop centrifuge at room temperature. Our favorite mini benchtop centrifuge reaches 7,000 rpm, but if yours is a different speed just keep spinning until the liquid passes through to the collection tube. It doesn't take much centrifugal force since there aren't any filters.
- 2) For the best results, keep the beads suspended during all incubation steps. We do this using a 360° rotisserie, but you can also use a carousel to keep the beads suspended in solution. We do not recommend shaking the beads, but if you do make sure you use **gentle** agitation.
- 3) The capture tube has two holes near the top that prevent loss of liquid when closing the tubes. Make sure you don't cover these holes while closing the cap!
- 4) Try to avoid touching the bottom of the ProMTag capture tubes, as anything that gets on the bottom of the tube may end up in your final sample.
- 5) We recommend briefly vortexing the ProMTag capture tubes at multiple points throughout the protocol to aid in resuspension of the beads. **Never vortex the capture tube alone**. Always vortex using the 2 mL waste tube as an adapter to avoid touching the bottom of the capture tube.
- 6) We do not recommend pipetting to mix at any stage where the capture resin is present, as the resin will stick to the tip and result in suboptimal yield. Try to avoid touching the resin with a tip at any point.

# Protocol for protein sample extraction, cleanup, and digestion using the UPECK-Pep kit

- 1)** If you are starting with a prepared lysate or proteins, skip this step but be sure to read the note above regarding lysis buffers. If you are starting with intact cells, start here. This is our favorite protocol to prepare a cell lysate, but if it won't work for your sample contact us and we'll help you figure out an alternative.
  - a.* Add LB to your pelleted cells. The volume you add is up to you, but we recommend adding about 2x the volume of your pellet (Ex: if your pellet is 100  $\mu$ L, add 200  $\mu$ L LB).
  - b.* Incubate the cells in a boiling water bath for 15 minutes.
  - c.* Optional: Sonicate (We use a Branson Probe sonicator) set at 30% power, 30% duty cycle for 10-15 pulses until the viscosity dropped for particularly stubborn cells.
  - d.* Pellet cell debris by spinning at  $\sim$ 3300 xg for 10 minutes.
  - e.* Transfer the supernatant to a fresh tube and assay protein concentration via your favorite method. We recommend a Bradford or BCA assay.
- 2)** Add 50  $\mu$ g to 200  $\mu$ g of protein in lysis buffer to a protein low bind tube provided with the kit. Bring concentration of your protein sample to 1 mg/mL by diluting using the provided LB.
  - a.* If your sample is dilute, no worries! You can use a larger starting volume if necessary. However, we recommend keeping the input volume under 200  $\mu$ L if possible.
- 3)** Add 10  $\mu$ L DTT and incubate your sample at 56°C for 30 minutes.
  - a.* The first time you use the kit, resuspend the provided DTT in 100  $\mu$ L water. After resuspending the DTT, store it at -20°C.
- 4)** While your sample is incubating, resuspend one tube of IAA in 20  $\mu$ L water. Vortex or pipet thoroughly to mix.
- 5)** Add 10  $\mu$ L IAA and incubate your sample at room temperature in the dark for 30 minutes.
- 6)** Add 17  $\mu$ L of ProMTag to your protein sample and pipette up and down or vortex briefly to mix well. Spin briefly to collect your sample at the bottom of the tube if necessary.

- 7)** Incubate for 30 minutes at room temperature to tag your proteins.
- 8)** During the last 5 minutes of the incubation, prepare the ProMTag capture tube for use.
  - a. Take a capture tube from the bag, remove the parafilm, and place it into one of the provided 2 mL waste collection tubes.
  - b. Spin briefly (~20 seconds) in a tabletop centrifuge until all the liquid has passed into the waste collection tube. Discard the flow through. Spin again if necessary to remove all liquid. The resin will be bright white when they're dry.
  - c. Add 300  $\mu$ L of WB3 buffer and vortex briefly (~2 seconds).
  - d. Spin briefly in a tabletop centrifuge until all the liquid has passed into the waste collection tube. Discard the flow through.
  - e. Reapply parafilm to the bottom of the tube
- 9)** Once step 7 and 8 are complete, bring the final volume of your tagged protein sample to 200  $\mu$ L with the provided lysis buffer.
- 10)** Add the tagged protein sample to the ProMTag capture tube. Vortex briefly to mix.
- 11)** Incubate the tube at room temperature with gentle rotation for 15 minutes. We use a rotisserie, but you can use anything that works to keep the resin suspended, so long as it is gentle (no harsh shaking). As the reaction proceeds, the pink solution should turn colorless.
  - a. If your sample is still pink after 15 minutes, allow it to incubate for a few more minutes.
- 12)** Remove parafilm spin briefly in a tabletop centrifuge until all the liquid has passed into the waste collection tube. Discard the flow through.
- 13)** Add 300  $\mu$ L WB2 to the sample. Vortex for ~5 seconds to mix. Spin briefly in a tabletop centrifuge until all the liquid has passed into the waste collection tube. Discard the flow through.
- 14)** Repeat step 13 one time.
- 15)** Add 300  $\mu$ L WB3 to the sample. Vortex for ~5 seconds to mix. Spin briefly in a tabletop centrifuge until all the liquid has passed into the waste collection tube. Discard the flow through.
  - a. Beads may stick along the sides of the capture tube after this wash. This is normal. Do not attempt to scrape the beads down with a pipet tip, etc.

- 16)** Add 300  $\mu\text{L}$  deionized water to the sample. Vortex for  $\sim 2$  seconds to mix. Spin briefly ( $\sim 10$  seconds) in a tabletop centrifuge until all the liquid has passed into the waste collection tube. Discard the flow through.
- 17)** Repeat step 16 one time.
- 18)** Reapply parafilm to the bottom of the tube
- 19)** Add 200  $\mu\text{L}$  EB to the capture resin. Briefly vortex the sample, then gently flick or tap the capture tube to bring most of the liquid back to the bottom of the tube.
  - a. From this stage until elution **do not** centrifuge the sample to bring the liquid back to the bottom of the tube, as doing so will result in loss of proteins or peptides.
- 20)** Incubate the tube at room temperature with gentle rotation for 15 minutes.
- 21)** Add 200  $\mu\text{L}$  mP-Trypsin to your sample and incubate at  $37^\circ\text{C}$  for 1 hour.
  - a. Longer digestion times will not improve yield. We recommend that digestion does not exceed 1 hour.
- 22)** Remove parafilm and transfer the ProMTag capture tube to a provided low bind protein tube.
- 23)** Spin briefly in a tabletop centrifuge until all the liquid has passed into the low bind protein tube. **DO NOT DISCARD THE FLOWTHROUGH.**
- 24)** Add 200  $\mu\text{L}$  EB to the capture resin. Vortex for  $\sim 5$  seconds to mix.
- 25)** Incubate the tube at room temperature with gentle rotation for 15 minutes.
- 26)** Return the capture tube to the same low bind protein tube. Spin briefly in a tabletop centrifuge until all the liquid has passed into the low bind protein tube.
- 27)** That tube now contains your sample of pure peptides in an acidic, volatile buffer. If desired, you can concentrate the sample by drying in a speed vac, or if you would like to use it as is, be sure to neutralize the sample with your buffer of choice.