

# perio care®

A highly elastic periodontal dressing which sets resiliently hard, does not chip, and provides durable protection of tissue.

## periodontal dressing

### Directions

**Helpful Hint:** To accelerate setting time, add 1-2 drops of water to PerioCare® paste and gel and incorporate water into the mix, or mix PerioCare paste and gel using a wet spatula.

1. **Extrude** equal amounts of paste (tube 1) and gel (tube 2) onto mixing pad.
2. **Mix** paste and gel together, using spatula, until color becomes uniform throughout. Thorough blending usually takes about 20 seconds.
3. **Form** PerioCare mixture into log shape on mixing pad, using spatula.
4. **Test** for tackiness with spatula or water-moistened finger. PerioCare is ready to use with moistened fingers in about 45-60 seconds. (you may notice that PerioCare is exothermic during this period.)
5. **Mold** PerioCare into desired shape using fingers. PerioCare will remain workable for approximately 4-5 minutes. If additional working time is required, place PerioCare in hot water to soften.
6. **Clean** spatula and instruments immediately with Wonder Orange™ orange solvent or with alcohol.

## Pulpdent perioCare Periodontal dressing

Rx only 

- Easy to measure: use equal amounts of paste and gel.
- Easy to mix: paste and gel combine easily and completely.
- Complete: no additional ingredients are necessary.
- Ready to use quickly: 45-60 seconds
- Readily formed into cords of any thickness or length for easy application in the mouth.
- Cohesive: PerioCare melds to itself during working time.
- Adhesive: PerioCare stays in the mouth as long as needed without shrinkage or disintegration.
- Easily removed.
- Naturally colored.
- Neutral aroma and taste.
- Contains no eugenol.
- Safe: contains no asbestos.
- The package is designed for efficient and convenient office storage. Each tube is clearly marked to avoid confusion between paste and gel.
- Economical for the dentist. Package contains 6 fl. oz. (177.4 ml). A mixing pad is included in the PerioCare kit.



### PULPDENT Corporation

80 Oakland Street • Watertown, MA 02472 • U.S.A.  
(617) 926-6666 or (800) 343-4342 • Fax (617) 926-6262  
www.pulpdent.com • pulpdent@pulpdent.com

XP-PC-IN-06w  
REV: 03/2014