

Bullseye Safety LANCETS™

Phoenix Healthcare Solutions is proud to introduce the Bullseye Safety Lancet™. The Bullseye is the first retractable lancet with a one-of-a-kind "Color Indicator" that disappears after the device has been used. The Bullseye Safety Lancet™ offers patients and health care professionals unrivaled safety and security.

Features & Benefits

Safety

- Color Indicator – Provides a visual to let you know if the Lancet has been used
- Single use retractable Lancet eliminates cross contamination




Comfort

- Ergonomically Designed
- Quick spring action improves patient comfort
- Tri-Beveled Needle Tip

Easy To Use

- Pre-cocked & Ready To Use
- Two Simple Steps: Remove cap & press button to fire



Bullseye Colour	Product Type	Puncture Depth	Needle Diameter	Blood Flow
	High-Flow Bullseye	2.2 mm	21G Needle	100 - 175ul (Approximately)
	Normal-Flow Bullseye	1.8 mm	26G Needle	40 - 60ul (Approximately)
	Low-Flow Bullseye	1.6 mm	28G Needle	20 - 40ul (Approximately)



+1-954-880-4272
mikem@phoenix-hs.com

2890 W State Road 84 | Suite 104, Fort Lauderdale, FL 33312 | USA
T + 1.954.607.6080 | F + 1.954.376.8203
www.phoenixhealthcaresolutions.com

