













Symbols Glossary

SYMBOL	REFERENCE No.	SYMBOL TITLE	EXPLANATORY TEXT
Symbols from <i>ISO 15223-1:2016 Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements</i>			
	5.4.2	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	5.2.3	Sterilized using ethylene oxide	Indicates a medical device that has been sterilized using ethylene oxide.
	5.1.5	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	5.1.6	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	5.1.3	Date of manufacture	Indicates the date when the medical device was manufactured.
	5.1.1	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directive 90/385/EEC, 93/42/EEC and 98/79/EC.
	5.1.4	Use-by-date	Indicates the date after which the medical device is not to be used.
	5.2.8	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	5.3.7	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.
	5.4.3	Consult instructions for use	Indicates the need for the user to consult the instructions for use.

Symbols Glossary

SYMBOL	REFERENCE No.	SYMBOL TITLE	EXPLANATORY TEXT
Symbols from <i>ISO 15223-1:2016 Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements</i>			
	5.5.1	In vitro diagnostic medical	Indicates a medical device that is intended to be used as an in vitro diagnostic medical device.
	5.5.5	Contains sufficient for <n> tests	Indicates the total number of IVD tests that can be performed with the IVD.
Symbols from <i>21 CFR 801.109</i>			
Rx Only	N/A	Prescription Use Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.