

## Comparison of the Displacement of the MEDmetric® KT 2000™ to the BBM Knee Arthrometer (KA)

Performed by Blue Bay Medical, Inc. (BBM) 2/26/2018 Using Gauge Blocks

Gauge Block		Knee Arthrometers	
in	mm	KT	KA
0.00	0	0.0	0.0
0.20	5.08	5.0	4.9
0.40	10.16	10.0	10.0
0.60	15.24	14.6	15.1
0.80	20.32	19.3	20.1
1.00	25.4	24.1	24.8

