


THE MAGIC TRICK THAT BEATS YOUR DRUG TEST



Nona Scientific vs. Standard Lab Tests (NonaClear 2.0)

The new science of specimen validity testing. Give us a call at **352.260.4700** or email **info@nonascientific.com** to trial NonaClear 2.0 specimen validity testing today.



Here's How They Are Tricking You

1. Present synthetic urine as a substitute for a clean sample.
2. Add synthetic urine to actual urine to cut samples to RX levels within prescription range.
3. Add "Spike" to synthetic urine to show positive for RX, while masking abuse or presence of other drugs.
4. "Reverse Flow" catheterization of synthetic urine into subject for a clean observed sample.

Test results of synthetic, adulterated, and authentic urine samples analyzed at the Forensic Toxicology Drug Testing Laboratory, Fort Meade, Maryland, USA

Sample No.	Description	Creatinine (mg/dL)	Specific Gravity	pH	Oxidant Activity (ug/mL nitrite equiv)	Interpretation
1	UPASS	92.6	1.013*	7.2	0	Valid
2	Synthetix5	21.4	1.008*	7.4	3	Valid
3	Dr. Greens	61.8	1.011	5.4	0	Valid
4	UTAK	25.9	Not tested	3.63	0	Adulterated
5	Liquid Gold	91.4	1.012*	6.5	0	Valid
6	Whizzinator	9.6	1.0014	8.1	0	Valid / Dilute
7	Xstream	58.9	1.009*	6.8	0	Valid
8	Liquid Gold	36.8	1.009*	7.8	0	Valid

*Measurements taken separately at Tripler Arm Medical Center Clinical Core Laboratory.