



**6" x 6" Gel Pillow**  
with white cotton cover

**6" x 9" Gel Pillow**  
with printed cover

# Dandle•LION Gel Pillows

With Disposable, White or Printed Cotton Covers

| Size                                         | Description                     |
|----------------------------------------------|---------------------------------|
| <i>(Standard thickness: ~5/8")</i>           |                                 |
| 3" x 6"                                      | Small Roll                      |
| 6" x 6"                                      | Square Pillow                   |
| 6" x 9"                                      | Also fits XS Dandle ROO or ROO2 |
| 6" x 11"                                     | Also fits S Dandle ROO or ROO2  |
| 6" x 12"                                     | Also fits M Dandle ROO or ROO2  |
| <i>(Pediatric or Adult thickness: ~3/4")</i> |                                 |
| 6" x 7"                                      | Pediatric Gel Pillow            |
| 6" x 10"                                     | Pediatric Gel Pillow            |
| <i>(Pediatric or Adult thickness: ~1")</i>   |                                 |
| 8" x 13"                                     | Pediatric Gel Pillow            |

DandleLION soft gel-filled pillows are designed to provide a comfortable cushion to meet the support needs of micro-preemie, preemie, and full-term infants. Gel pillows can help reduce the need for frequent re-positioning and make immobilized babies more comfortable. They can also be used as a mattress for very low birthweight babies or as a prone roll for bigger babies.

- Three covers available: white cotton, printed cotton, and disposable
- Non-toxic, water based food gel
- Guaranteed for at least three months
- Guaranteed mold resistant
- DEHP and BPA-free, and not made with natural latex rubber
- Single patient use (unless hospital validates their own cleaning process)

**Ordering Information**  
See Electronic Catalog on resource page of website

Made in the USA

Dandle•LION Medical | 22 Shelter Rock Lane | Danbury, CT 06810 | 203-791-9001 | www.dandleLIONmedical.com

Rev5:0620

