

Post Mastectomy Supply Use & Care Instructions

Proper care/maintenance/disinfection of mastectomy bras

Hand wash in warm water with a mild soap. No bleach. Line dry.

Proper care/maintenance/disinfection of leisure forms

The leisure form is intended for temporary or occasional use only and is to always be used in conjunction with a proper cover. Do not wear this product directly on an open wound. To clean, hand wash foam portion and protective cotton cover separately in warm water after each use. Allow to air dry. Avoid excessive heat. Do not machine wash or tumble dry.

The materials contained in this product may be sensitive to x-ray equipment. When traveling, please pack in your checked-in luggage. This product may contain lead. Please wash hands after handling if the cover is damaged and do not allow small children to handle the leisure form. Please dispose of your breast form in accordance with local and national laws and regulations.

Proper care/maintenance/disinfection of silicone forms

Treat your breast form as you would your own skin/body. Take care with sharp objects. Finger nails, brooches, scissors and pets claws can cause permanent damage to the breast form's film. Do not pinch, rub or wring the breast form at any time. Avoid applying powder, perfume, deodorant, cream, moisturizer, tanning products or any lotions to the breast form, as well as to the chest wall. Please use care when inserting the breast form into the protective cover as repeated, unnecessary squeezing or pulling of the breast form can cause damage.

Wash breast form by hand daily with a mild soap and lukewarm water, rinse thoroughly, and then towel dry. Do not wring. Store the breast form in its original packaging and away from direct heat when not being worn. This allows the breast form to maintain its shape and protect it from accidental harm.

In the unlikely event the breast form contains small air bubbles, please allow a few days for them to dissipate. Air bubbles do not affect the integrity or performance of the breast form.

A breast form should only be worn when the surgical site is completely healed and radiation treatments are finished. Please consult with your physician before wearing a weighted breast form. Do not wear if skin is broken or irritated. If any irritation occurs, discontinue use immediately and contact your physician.

There is no guarantee that the adhesive properties of an attachable prostheses will work equally for all women. The durability of the gel adhesive is directly dependent on the level of care by the wearer and their commitment to clean and maintain the gel. The adhesive is excluded from the two year warranty. There is a limited two year warranty on breast forms. Replacement is for manufacturer's defects only, not consumer abuse, neglect or loss. Failure to clean breast forms will void your warrant. Replacement breast forms are from date of original purchase