

Using a Tablet to Sign Documents and Forms in Tracker

This document outlines the setup and process for signing documents and clinical forms using a Windows Surface Pro tablet.

Requirements

Unlike with our previous methods of incorporating signatures, when you are using the Windows Surface Pro tablet there is no additional hardware required. Instead the *office* will be required to purchase the Topaz drivers directly from Topaz and install them on the appropriate machine

NOTE: This driver is NOT provided by The Bridge Network and is licensed **directly from Topaz.**

Offices that want to use this must [contact Topaz and request the SigPlus Pro](#) for a Windows 8 Professional Tablet.

1. SigPlus Pro for Tablets
2. Tracker 11.18 or higher
3. MS Office plugins and / or Adobe Plug Ins as needed

Once the SigPlus Pro for Tablets is installed on your system, the stylus can immediately be used in the clinical forms. Additional plug-ins can be used for signing in additional applications.

Troubleshooting Tips

The DemoOCX program is also installed with the Tablet drivers. This program can be used to test the drivers and ensure that the machine can read / pick up signatures in the system.

Note that you CANNOT install the Tablet PC driver and the USB / Bluetooth driver on the same machine as it uses the same ini file. The driver that was installed last will be the one that machine is using.