



www.koomedical.com

CPAP Masks

PEEP valvePEEP

- Adjustable from 5 to 20 cmH₂O for CPAP ventilation.
- With connection OD 33 mm for CPAP masks KM-301A, KM-301DA, KM-303
- With two reducer connectors: 1x22Mx30,5 mm ; 1x22Mx33mm they allow a universal use of the valve for helmets, masks and circuits.
- Single patient use.
- Code: KM-809 (50 items per box).
- Each valve is packed individually with a description of the contents.



Code	Description	Dia. Connector	Weight g	Items per Box Pcs
KM-809	valvePEEP	33mm OD	22	50
KM-809A	valvePEEP	22mm OD	22	50

PEEP stepPEEP



- 7 discrete steps from 5 to 20 cmH₂O for CPAP ventilation
- With connection 33 OD mm for CPAP masks
- KM-301A, KM-30DA, KM-303A.
- KM-810 includes two reducer connectors : 1x 22 M x 30,5 mm ;
- 1x 22 M x 33 mm allow for universal
- use of the valve for helmets, masks and circuits
- Single patient use

Code	Description	Dia. Connector	Weight g	Items per Box Pcs
KM-810	stepPEEP	33mm OD	22	50
KM-810A	stepPEEP	22mm OD	22	50

KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for EMEA (Europa Middle East Africa) Markets - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com
Responsible for Asian Pacific Market - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com