



www.koomedical.com

Respiratory Care

Tri-Gym & Flo-Gym

Lung Exercisers

- In polycarbonate, tube and mouthpiece in polyethylene
- Latex free
- For inspiration and expiration
- With three measuring cylinders for the control of the flow
- Flow control is done through the three colored balls, which indicate the amount of flow per second 600 cc green, blue 900 cc, 1200 cc red.
- With corrugated tube and mouthpiece
- Single patient use



Tri-Gym

Code KM-805

a single packaging, 24 items per box.

Flo-Gym

Can be disassembled, cleaned with mild detergent and reassembled by hand.

Code: KM-802



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for EMEA (Europa Middle East Africa) Markets - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com
Responsible for Asian Pacific Market - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com