

# NIV masks

Non-Invasive Ventilation Face mask

## Bluestar

*For lung ventilation devices. Modality CPAP and BILEVEL assisted and/or supervised.*



- Inflatable cushion adaptable to the patient's face with a firm and comfortable fit. No skin breakdown.

### The kit is provided with:

Transparent PVC mask with inflatable cushion, valve with 6 grips to fix the mask, adapter connector (KM-1019) with ports for O<sub>2</sub> and CO<sub>2</sub>, pressure monitoring and a double caps for the introduction of a nasal cannula, DEHP and BPA free.

- Pump to inflate and deflate the cushion.
- Headgear's cuff with a rapid removal system.
- Latex Free.
- Possible to link with connector (KM-903 not included) provided with a hole for the CO<sub>2</sub> discharge and a port for the improvement of O<sub>2</sub>.
- Single patient use.



*The kit is inside a box with a label that indicates the contents.*

Code	Description	Size	Weight g	Dead space ml	Dia. of the connector	Items per box Pcs
KPM-001L	face mask NIV	extra large	91.0	290	22F	50
KPM-001M	face mask NIV	large	79.5	240	22F	50
KPM-002M	face mask NIV	medium	70.0	220	22F	50
KPM-003M	face mask NIV	small	65.0	175	22F	50
KPM-004M	face mask NIV	pediatric	46.0	85	22F	50
KM-1019	adapter connector	standard	19.0	17	22M	100

### KOO MEDICAL EQUIPMENT

Web: [www.koomedical.com](http://www.koomedical.com)

Responsible for EMEA (Europa Middle East Africa) Markets - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: [info@kooeurope.com](mailto:info@kooeurope.com)  
Responsible for Asian Pacific Market - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: [sales@koo-asia.com](mailto:sales@koo-asia.com)