



Proven Nipple Aspirate Collection System

Helps Identify Women at High Risk for Breast Cancer.




*F*act is all women are at risk of developing breast cancer. detection of molecular changes within breast ducts helps you and your patients take control of breast health.

Now you can easily identify which of your patients are at the greatest risk of developing breast cancer using a simple and reliable breast fluid test in your office, and it takes only minutes.

The HALO™ Nipple Aspirate Fluid Collection system is the first and only fully automated, noninvasive method for collecting ductal fluid for proteomic analysis. It's an easy and efficient way to routinely evaluate and manage your patients' breast health.

ABNORMAL CELL GROWTH DETECTION



*If breast cancer is discovered early, while contained within the breast ducts (Stage 0 or 1), survival rate is almost 100%. Earlier detection could ultimately reduce breast cancer mortality.

*reference:
American Cancer Society
www.cancer.org/docroot/CRI/content/CRI_2_2_3X_After_the_tests_Staging_5.asp?site=area=
 Accessed 03/13/07





System Components



HALO Console

- Uses a combination of warmth, suction, and massage to elicit Nipple Aspirate Fluid (NAF)
- Noninvasive, 5-minute cycle is fully automated and requires no specialized training or certifications
- Microprocessor-controlled, self-monitors to assure every cycle remains within specifications for vacuum, temperature and compression

Not indicated for use in lactating women.

The HALO system has not been tested on pregnant women.

These products and their use are covered by one or more of U.S. Patent Nos. 6,670,610, 6,712,785, 6,899,696, 6,866,994, 6,517,513, 6,875,184, 6,981,950.



HALO Patient Interface Cassette

- Adjustable Breast Cups and fluid reservoir with connecting tubing to conduct suction, warmth, and massage to breasts
- Adjustable Breast Cups fit virtually any breast size
- Replace every 3 to 4 months
- Simple to install



HALO Disposable Sample Collection Cups

- Soft, flexible, latex-free
- Prevents patient contact with HALO System so equipment cleaning is not required
- Includes 2 sample collection swabs and prep pad



HALO Cart (optional)

- Holds HALO console securely for room to room movement
- Locking wheels

Catalog Number	Description	Quantity
NAFD-100	HALO Console	1 each
NAFCART-01	HALO Cart	1 each
NAFCASS-01	HALO Patient Interface Cassette	1 each
NAFCUPS-01	Disposable Sample Collection Cups	24 pairs/case

